


Antepartum Fetal Monitor BT-330

BT330 is optimized Color LCD fetal monitor to check your fetal health

- 7" Wide TFT-Color LCD Display
- Up to 150 patients data saving
- Wide range of FHR : 30~240bpm
- Simultaneous Fetal Heart sound Play and Record
- Easy User Interface



Fetal Monitor **BT-330**

Specifications

FHR	<ul style="list-style-type: none"> • Dual Pulsed Doppler Ultrasound • Intensity : 10mW/cm² or less • Accuracy : $\pm 2\%$ of Range 	<ul style="list-style-type: none"> • Frequency : 0.984 MHz • FHR Range : 30~240bpm FHR
FM	<ul style="list-style-type: none"> • Auto-detection of Dual Fetal Movement 	
UC	<ul style="list-style-type: none"> • External Type • Frequency Response : DC~0.5Hz Reference (zero) Control 	<ul style="list-style-type: none"> • Measurement Range : 0~99 units
Printer	<ul style="list-style-type: none"> • Thermal Array Type • Auto Record Period : Off, 10, 20, 30, 40, 50, 60min 	<ul style="list-style-type: none"> • Print Speed : 1,2,3 cm/min, high speed • Paper Feeding Function
Display	<ul style="list-style-type: none"> • 7" Wide TFT-Color LCD : 800(W) x 480(H) 	<ul style="list-style-type: none"> • 3Channels (FHR I, FHR II, UC)
Sound	<ul style="list-style-type: none"> • Doppler Sound Alarm 	<ul style="list-style-type: none"> • Sound Information Sound
Set-up	<ul style="list-style-type: none"> • Alarm Upper / Lower Value, Time & Date Print Speed • Trend [Data Saving for 450 hours(3hr/Person)], Name ID 	
Function	<ul style="list-style-type: none"> • Mark Function • Auto Print Function • Heart Sound Recording in BCM330 	<ul style="list-style-type: none"> • Display/Printer Zoom-in Function • FHR II Offset Function(DOP +20bpm)
Power	<ul style="list-style-type: none"> • Input : AC 100~240 V (50/60 Hz) • Battery Operation : continuously 2 hours 	<ul style="list-style-type: none"> • Output : DC 16V, 2.8A Consumption
Dimension	<ul style="list-style-type: none"> • 191mm(H) x 187mm(L) x 201mm(D) 	
External Link	<ul style="list-style-type: none"> • Wire : RS-232C 	<ul style="list-style-type: none"> • Wireless : Bluetooth
Standard Accessory	<ul style="list-style-type: none"> • Ultrasound Doppler Probe 2ea • Event Marker 1ea • Power Cord 1ea • Operation Manual 1ea 	<ul style="list-style-type: none"> • UC Probe 1ea • Print Paper 2ea • Ultrasound Gel 1ea • Probe Belt 3ea • AST(Acoustic Stimulator) 1ea(Optional)
Warranty	<ul style="list-style-type: none"> • 2 Years(Probe excluded) 	

PC Interface

BCM-330

(Central Monitoring System)

- Remote control up to 8 monitors by central PC
- Simultaneous Fetal heart sound Play and Record
- Unlimited Multi language support
- Wire(RS-232C) and Wireless(Bluetooth) connection

